

IO Endo | rotary file



Compatible with
all endodontic
micromotors
on the market

AND IF YOU WERE TO CHANGE TO FILES OFFERING REAL PROGRESS IN TERMS OF CONVENIENCE AND FLEXIBILITY?

With its **IO Endo files**, IONIX is offering an innovative product that adapts perfectly to all types of canal: simple or complex, very narrow or curved

performance

IO Endo files adapt perfectly to the shape of the canal allowing you to work safely. The fineness of its non-cutting tip combined with the two cutting edges of its blade ensure excellent cutting efficiency, thereby avoiding any risk of obstruction, deviation or perforation.

efficiency

The IO Endo range offers you a complete solution: IO Endo Small, Medium or Large, for the simplest to most complex root canal treatments; IO Endo Retreatment files for removing filling material.

simplicity


Each pack of IO Endo files contains a trio or a quartet of files that are very easy to identify thanks to the coloured rings etched onto the shaft. Save time when carrying out your dental procedures. Using IO Endo files ensures you can work quickly and efficiently, and most importantly, adapt to each specific case.

YOUR PARTNER FOR OVER 25 YEARS



Founded in 1993 and established in France in the Bordeaux region, the company, **IONIX designs and manufactures specialised endodontic devices**. Intended for general practitioners devoted to improving the quality of their treatments as well as the comfort provided, **the entire IONIX was designed with everyday practice and ease of use in mind**.


With a worldwide presence through its powerful network of distributors, **IONIX is one of the most widely recognised specialists in the dental sector**.



IO Endo


rotary file

TECHNICAL CHARACTERISTICS



IO Endo Small
(narrow apical)

For very narrow or curved canals, to be used in addition to initial preparation with IO Endo Medium.



IO Endo Medium
(normal apical)

For straight or slightly curved canals with a standard diameter.



IO Endo Large
(wide apical)

For wide canals, to be used in addition to initial preparation with IO Endo Medium.



IO Endo Retreatment

For removing old filling material.

DUE FOR
RELEASE
IN 2019

ENDY 7. THE NEW GENERATION OF MICROMOTORS WITH INTEGRATED APEX LOCATOR

Speed from 200 to 1000 rpm - Torque from 0.5 to 3.5 N.cm



All our IO ENDO files are compatible with the endodontic micromotors currently on the market.

The Endy 7 is an upgrade of the Endy 6200. Scheduled to be launched in 2019.

A REAL PROGRESS IN ENDODONTIC PRACTICE

ionyx

34 rue de Fleurenne
F-33290 BLANQUEFORT

T. +33 5 56 35 34 27
F. +33 5 56 95 16 20

ionyx@ionyx.eu
www.ionyx.eu

